

Upgrade to new benefits

Introducing the CORE 2 Console

Today's next-gen console combines Stryker's innovation with your real-world insights, providing you with more efficiency in your small bone cases.



	Original CORE Console	CORE 2 Console
Boot up time		2x faster
Screen view	6.5" with 107° approximate viewing angle	7" with 170° wide view angle
Color	18 bit	24 bit with graphics
Larger icons		✓
Flatter menus/page structure		✓
Illumination rings show matched hand pieces/foot pedals		✓
USB transfer of user preferences		✓
Methods to adjust settings	2	3 – new vertical touch screen scroll bar
Connector port design	Recessed	Flush with surface
Connector port location		Flipped; more space around irrigation port
Audio alerts		More; with distinct tones for each type
Error codes	Numeric code	Numeric, written instruction and rep contact stored internally
Servicing and stored settings	User presets lost; re-entry required	User presets saved; no re-entry needed
Dimensions	12.5" w x 5" h x 17" d	13" w x 5.4" h x 17.4" d
Weight	9.07 kg	7.8 kg
Languages	5	18

Upgrading is easy

Your Surgical sales representative is here to make your upgrade quick and easy so you'll get the most from your system right out of the gate. For more details, contact your Surgical sales representative or **visit strykermeded.com**.

CORE 2 Console

5400-052-000

Accessories:

Single handpiece irrigation cassette (disposable)

5400-050-001

Dual handpiece irrigation cassette (disposable)

5400-050-002

Handpiece cord

5100-004-000

Accessory Cart 5400-410-000

Bi-directional footswitch

5100-008-000

NSE footswitch 5400-007-000

RemB Handpieces

RemB Micro Drill

6400-015-000

RemB Oscillating Saw

6400-031-000

RemB Sagittal Saw

6400-034-000

RemB Recipricating Saw

6400-037-000

RemB Universal Driver

6400-099-000

RemB Dedicated Wire Driver

6400-062-000

RemB Handswitch

6400-009-000

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: CORE, RemB, Stryker. All other trademarks are trademarks of their respective owners or holders.

9100-004-827SSP

G55/PS

Copyright © 2017 Stryker

Printed in USA

Stryker Australia Pty Ltd
8 Herbert Street,
St Leonards NSW 2065
T: 61 2 9467 1000
F: 61 2 9467 1010

Stryker New Zealand Limited
515 Mt Wellington Highway,
Mt Wellington Auckland
T: 64 9 573 1890
F: 64 9 573 1891

www.strykermeded.com